


# Long Life Deuterium UV Detector Lamps

Over 30% cost saving with no compromise on OEM Quality



- Premium Quality Long life Lamps
- 2000 hours lifetime guarantee
- OEM compatible lamps from worlds leading manufacturers
- Every lamp Individuality tested and certified
- Lamps available for all leading makes of detectors



## Premium Quality- Available Stock- OEM Level Technical Support

Neta works in the manufacturers/ pre-aligns of Deuterium lamps for independent instrument service companies. In recognizing the importance of using lamps supplied by the OEM's, Neta has formed strong partnerships with manufacturers who produce the majority of the lamps fitted in today's leading detectors, including Shimadzu and Waters models. Bare lamps supplied directly to manufacturer are precision aligned and tested y in accordance with stringent ISO 9000 certified procedures.

All Neta lamps are certified recording the date of manufacture, serial number, tester's name and the intensity counts reading for future traceability. Only lamps made to the strictest tolerances receive the Neta Guaranteed Lamp Seal. This means that Neta provides lamps without any compromise on quality compared with OEM specifications. We offer a comprehensive range the Shimadzu SPD10A and Waters 2487. A range of lamps are kept in stock in the factory and our sales offices, close to our customer ensuring that when replacements are needed there is no delay.

To assist our customers with pre-sale information and after sales-support Neta is able to offer a new, higher level of technical support to our customers which was previously only available to the instrument OEMs.

Premium quality lamps, available from stock with the highest level of technical support before and after the sale.



**Shimadzu SPD10A Long Life D2**



**Gilson 115/116/117/118/119/151/  
152/155/156**



**Gynkotek UVD  
320/160/170S/340 DAD**



**Waters 996, 2996 PDA Long Life  
Deuterium**



**TSP UV100/1000/2000/3000  
Linear Focus and 200 Series**



**Varian UV50/100/200 9050  
Prostar 310 D2 long life**



**Waters Alliance 2487 Long Life D2**

With our range of replacement Long Life Deuterium UV Detector Lamps we are able to supply the highest quality products available whilst offering great savings over lamps from the original equipment manufacturers (OEMs).

Our range of Long Life Lamps have a guaranteed life time of 2000 hours or 1 year, so not only do you get a great price, you also receive a comparable warranty to that offered by the OEM's providing complete piece of mind.

To ensure our quality we only take bare lamps from the same renowned manufacturers that supply the OEMs. *These lamps are then precision aligned and tested in accordance with stringent ISO 9000 accredited procedures. All lamps are certified, recording the date of manufacture, serial number and intensity count for future traceability.* Only lamps made to the strictest tolerances receive our Guaranteed Lamp Seal ensuring no compromise on quality compared to OEM specifications.

**\*\*\*Our Most Popular Long Life Lamps\*\*\***

<u>Part Number</u>	<u>Description</u>	<u>OEM List price</u>	<u>Net Sales Price</u>
NS-DIO-102LL	Dionex PDA-100 AD-25 D2 Lamp	>\$700.00	<b>\$500.00</b>
NS-GIL-101	Gilson 115/116/117/118/119/151/152/155/156 D2 Lamp	>\$450.00	<b>\$365.00</b>
NS-GIL-105	Gilson 170 D.A.D D2 Lamp	>\$650.00	<b>\$580.00</b>
LD-AGI-104LL	Agilent 1100 VWD Longlife D2 Lamp	\$630.00	<b>\$455.00</b>
LD-GYN-100LL	Dionex/Gynkotek UVD 320/340/160/170S DAD D2 Lamp	>\$550.00	<b>\$468.75</b>
NS-AGI-105LL	Agilent 1100, 1200 DAD Longlife D2 Lamp	\$920.00	<b>\$500.00</b>
NS-PER-100	Perkin Elmer Series 200 D2 Lamp	\$435.00	<b>\$335.00</b>
NS-PER-102	Perkin Elmer Series 200 DAD D2 Lamp	\$435.00	<b>\$430.00</b>
NS-SHI-101LL	Shimadzu SPD10A Long life D2 Lamp	\$705.00	<b>\$435.00</b>
NS-SHI-103LL	Shimadzu LC2010 Long life D2 Lamp	\$505.00	<b>\$390.00</b>
LD-TSP-102	TSP UV100/1000/2000/3000 Focus Spectrachrom D2 Lamp	>\$400.00	<b>\$325.00</b>
LD-VAR-103LL	Varian UV50/100/200 9050 Prostar 310 D2 Long life Lamp	\$715.00	<b>\$485.00</b>
NS-VAR-104	Varian Prostar 330 DAD D2 Lamp	\$715.00	<b>\$485.00</b>
NS-VAR-105LL	Varian Prostar 325 335 Long life D2 Lamp	\$750.00	<b>\$495.00</b>
NS WAT 102	Waters lamp 484/486	\$550.00	<b>\$435.00</b>
LD-WAT-104LL	Waters 996 2996 PDA Long life D2 Lamp	\$625.00	<b>\$435.00</b>
LD-WAT-105LL	Waters 2487 Long life D2 Lamp	\$700.00	<b>\$550.00</b>

**Telephone: (800) 343-6015      Fax: (609) 265-8213      Email: [rejean@netascientific.com](mailto:rejean@netascientific.com)**

**Neta is your single source for fast, accurate delivery  
of your high quality products anywhere in the USA.**

**For more information about our products**

- Visit us online at [www.netascientific.com](http://www.netascientific.com)
- Call: (800) 343-6015, Fax: (609) 265-8213
- Email: [rejean@netascientific.com](mailto:rejean@netascientific.com) or  
[orders@netascientific.com](mailto:orders@netascientific.com)



Notes:

---

---

---