



# FREE BEADS!

Purchase 5 packs of prefilled bead homogenizer tubes, get one FREE!

Benchmark Scientific offers a wide range of bead sizes and materials to meet your sample needs. Our Triple-Pure™ line is certified nuclease and protease free and tested to be free of stray DNA.

The prefilled, 2.0ml screw capped tubes are supplied in packs of 50 and can be used with the Benchmark line of homogenizers as well as most other bead homogenizers on the market.

## Prefilled Bead Tubes

	Catalog no.	Size/Material	Application	List Price
Triple-Pure™ Zirconium	D1032-01	0.1mm Zirconium	Bacterial lysis (gram +/-), small yeasts	\$146.00
	D1032-05	0.5mm Zirconium	Yeast, algae, spores	\$146.00
	D1032-0105	0.1/0.5 Mixed Zirconium	Bacterial lysis, yeast, algae, spores	\$146.00
	D1032-10	1.0mm Zirconium	Larger yeasts, algae, fungi	\$146.00
	D1032-15	1.5mm Zirconium	Soft tissues, liver, brain, adipose, spleen	\$146.00
	D1032-30	3.0mm Zirconium	Tough tissues, heart, muscle, leaves	\$146.00
	D1032-SK	Starter kit, 10 of each above	See above	\$146.00
Acid Washed/Heat Treated	D1032-60	1x6mm Zirconium Satellite	Dry grinding hard samples, seeds, bone, hair	\$146.00
	D1032-RF-60	As above, reinforced tube	Dry grinding hard samples, seeds, bone, hair	\$153.00
	D1031-01	0.1mm Silica/Glass	Bacterial lysis (gram +/-)	\$104.50
	D1031-05	0.5mm Silica/Glass	Yeast, algae, spores	\$104.50
	D1031-10	1.0mm Silica/Glass	Larger yeasts, algae, fungi	\$104.50
	D1033-28	2.8mm Stainless Steel	Hard samples, insects, tough plants	\$146.00
	D1033-30G	Garnet shards, 1x6mm zirconium	Skin, highly fibrous samples	\$146.00
	D1034-MX	0.1, 4mm glass, 1.5mm zirconium	Feces, environmental samples	\$146.00
	D1034-28	2.8mm stainless steel, 5ml tube	Larger, hard samples	\$195.00

All beads available in loose bulk also. Empty 2.0ml tubes, standard and reinforced, available. Call for details.

Free offer expires March 31, 2018. Use offer code FREEBEADS18. Valid on prefilled tubes only and all must be purchased at the same time. Can be mixed or matched.



### Benchmark Microtube Homogenizer Family

BeadBlaster™ 24

BeadBug™ 6

BeadBug™ 3

D2400	BeadBlaster 24, 24x2ml (or 12x5ml*)	\$6495.00
D1036	BeadBug™ 6, 6x2ml (or 2x5ml*)	\$1995.00
D1030	BeadBug™ 3, 3x2ml	\$749.00

\*with optional adapters