

# High-Speed

# DIAMOND MAX™ Centrifuge Tubes

**NEW!** Diamond Max™ Centrifuge Tubes are engineered for high-speed centrifugation and meet the demands of the most critical scientific work.

**BUY 2 GET 1 FREE!\***

### FEATURES AND BENEFITS:

**Ultra-High Performance** – High spin rates with maximum RCF ratings of up to 20,000 x g

**Premium Materials** – USP Class VI resin is also FDA CFR 21 compliant for food safety

**Certified** – Tubes are manufactured in a class 100,000 clean room, gamma irradiated to SAL 10<sup>-6</sup>, and are certified to be free from DNase and RNase, Human DNA, Pyrogens, and Endotoxins

**Leak Resistant** – Easy Grip Caps have an integral seal ring design for secure closures that eliminate seepage when used as recommended

**Value** – Efficient automated production reduces cost and assures high quality

**Ease of Use** – Exceptionally clear tubes have black printed graduation marks to full volume of 15mL or 50mL, a bright white marking space, and are packaged in heavy duty, zip, re-closable bags.

**Eco Friendly** – Paperboard racks are readily recyclable in the well-established paper recycling stream

*\*Buy any 2 cases of Diamond Max Centrifuge Tubes and get 1 case free. Offer valid April 1st 2018 thru June 30th 2018. Mix and match. Product of equal or lesser value free.*

**SPIN RATES  
UP TO 20,000 X G!**



Item #	Description	Max RCF (g)	Unit
6295	15mL Centrifuge Tube w/Cap, STERILE, 25/Bag, 20 Bags/Case	17,000	500
6296	15mL Centrifuge Tube w/Cap, STERILE, 25/Rack, 20 Racks/Case	17,000	500
6297	50mL Centrifuge Tube w/Cap, STERILE, 25/Bag, 20 Bags/Case	20,000	500
6298	50mL Centrifuge Tube w/Cap, STERILE, 25/Rack, 20 Racks/Case	20,000	500
6299	50mL Centrifuge Tube w/Cap, STERILE, Self-Standing, 25/Bag, 20 Bags/Case	7,000	500



**Neta Scientific®**

LAB SUPPLIES AND SOLUTIONS

[www.netascientific.com](http://www.netascientific.com)

4206 Sylon Blvd • Hainesport, NJ 08036

tel: 800.343.6015 • fax: 609.265.8213

# OFFICIAL REDEMPTION FORM



## Buy any 2 cases of Globe Scientific Diamond Max™ Centrifuge Tubes and get 1 case FREE!

**Please follow the instructions below to redeem promotions.**

1. Order products between April 1 and June 30, 2018.
2. Complete the form below, then e-mail this redemption form and proof-of-purchase (Neta Scientific invoice or packing slip) to Dennis Hoeft at [d.hoeft@globescientific.com](mailto:d.hoeft@globescientific.com)

Offers subject to change without notice. Please allow 4-6 weeks for delivery of free products. All redemption forms must be received no later than 30 days after the promotion's expiration date.

### Enter Account Information

Name \_\_\_\_\_ Date \_\_\_\_\_

Organization \_\_\_\_\_ Account No. \_\_\_\_\_

Ship to  
Address \_\_\_\_\_

Room/Bldg./Dept. \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

### Enter Offer Information

Item No. of	Qty.	Item No. of	Qty.

Attach additional sheet if necessary

Neta Scientific Proof of Purchase (invoice or packing slip) MUST accompany this redemption form!