



# NEW! **Locking** Microcentrifuge Tubes

**Available In 0.5mL, 1.5mL, 2.0mL & 5mL MACRO Tube**

New Globe Scientific Locking Microcentrifuge Tubes are engineered with a secure locking system that provides the security of a screw cap microtube with the simplicity of a snap cap tube.



*Tubes shown actual size*

## Features at a Glance:

**Secure** – Maintains sample integrity with locking cap which prevents accidental cap opening and evaporation during centrifugation

**Versatile** – Single use tubes for laboratory procedures including centrifugation, boiling, incubation, freezing and sample storage

**Quality** – Tubes are made from optimally clear USP Class VI medical grade polypropylene

**Certified** – Manufactured in a class 100,000 clean room and are certified to be free from RNase and DNase, Human DNA, Pyrogens, PCR Inhibitors, and ATP

**Ease of Use** – Tubes feature graduations, writing area, pierceable cap, are autoclavable, and are packaged in zip, re-closable bags

Item #	Description	Unit
111552LK	0.5mL, polypropylene (PP), natural,	500
111562LK	1.5mL, polypropylene (PP), natural	500
111572LK	2.0mL, polypropylene (PP), natural	500
111582LK	5.0mL, polypropylene (PP), natural	100



**Neta Scientific**  
LAB SUPPLIES AND SOLUTIONS

www.netascientific.com  
4206 Sylon Blvd  
Hainesport, NJ - 08036  
t 800.343.6015 | f 609.265.8213