



# Quality Biological<sup>®</sup> Endotoxin-Free Water

Eliminate the risk of toxic activity produced by bacteria that can interfere with cell structure and function.

<u>UOM</u>	<u>CAT. NO.</u>
500 mL	QB-118-325-101
1000 mL	QB-118-325-131

Click the part numbers to shop now



Order Now at 800-343-6015 or [www.netascientific.com](http://www.netascientific.com)

# Endotoxin-Free Water

---

Ideal for preparing solutions and reagents for transfection, transformation, cell culture and tissue culture.

Use when producing ultra-low endotoxin buffers and serum-free media.

Suitable when testing under pyrogen-free conditions (LAL testing) is required.

Reliable for rinsing and washing equipment, glassware and plasticware in pyrogen-free conditions.

- Mycoplasma tested
- DNase, RNase and protease tested
- 0.1 µm sterile filtered
- Meets < 0.005 EU/mL endotoxin by USP/BET-compliant Endosafe® nexgen-PTS assay
- Manufactured in an ISO 9001:2015 certified facility
- Storage conditions: 2°C-30°C

*This product is for research purposes only. Not for diagnostic use.*



**Neta Scientific®**  
LAB SUPPLIES AND SOLUTIONS



800-343-6015



sales@netascientific.com



4206 Sylon Blvd • Hainesport, NJ 20879



www.netascientific.com